Reserve double stage central unit

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These units are realised according to the **ISO 7396-1** and are suitable to be installed in certified medical gas pipeline systems.

The third source unit, according to the current regulatory standards, it must be able to supply gas to the pipeline system in case of malfunction of the primary sources.

The technical features of the reducers, valves, pressure gauges and switches and in general of all the components reflect the ones of the main unit as follow:

- First stage pressure reducer according to the **ISO 10524-2**;
- Shut-off valves with metallic seat and shutter; they have successfully withstood the adiabatic compression test as required by the **ISO 7396-1**;
- High and low pressure gauges according to **EN 837-1**;
- High pressure switches preset at 20bar, they have successfully withstood the adiabatic compression test;
- Gas specific inlet connections;

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Reserve double stage central unit **C-2210M-2250M-R** high flow rate

Technical Data

SUITABLE FOR THE FOLLOWING GASES: O₂, N₂O, Air, N₂, CO₂

MAX FLOW RATE: see flow rate diagram

MAX INLET PRESSURE: 200 bar (for O₂, Air and N₂); 110 bar (for N₂O and CO₂)

MIN INLET PRESSURE: 10 bar

OUTLET PRESSURE: 6 ÷ 10 bar (presetted at 6.5 bar)

STORAGE TEMPERATURE: -20 ÷ +60 °C

OPERATING TEMPERATURE: 10 ÷ 40 °C

OUTLET CONNECTION: soldering coupler Ø 28mm

PRESSURE RELIEF VALVE OUTLET CONNECTION: soldering coupler Ø 16mm

DIMENSIONS (mm): 330x230x210

FRAME MATERIAL: inox sheet

STANDARDS: ISO 10524-2; ISO 21969; ISO 15001; ISO 7396-1


ALL THE COMPONENTS OF THE SYSTEM SUBJECTED TO CYLINDER PRESSURE OR TO SINGLE FAULT CONDITION HAVE SUCCESSFULLY WITHSTOOD THE ADIABATIC COMPRESSION TEST IN ACCORDANCE WITH ISO 10524-2 AND ISO 15001

DESGREASED IN ACCORDANCE WITH ISO 15001

Gas Table	GAS	GAS SYMBOL	COLOR	HIGH PRESSURE CONNECTION
	Oxygen	O ₂	White	UNI 11144 (W 21,7 M x 1/14")
	Medical Air	AIR	White / Black	UNI 11144 (W 30 M x 1/14")
	Nitrous Oxide	N ₂ O	Blue	UNI 11144 (G 3/8" M)
	Carbon Dioxide	CO ₂	Grey	ISO 5145 (W 27 x 2 M)
	Nitrogen	N ₂	Black	UNI 11144 (W 21,7 F x 1/14")

Flow rate diagram

Pressure regulator 2250M | Test Gas N₂

Test gas Nitrogen
Inlet pressure 20 bar
Test temperature 15 °C
Conversion factors To find the flow of a different gas from the test gas, multiply the value read on the curve for the coefficient shown in the following table:

Air	0.941
Oxygen	0.935
Carbon dioxide	0.761
Nitrous oxide	0.795



Reserve double stage central unit **C-2210M-2250M-R** high flow rate

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