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AGSS
Station Back view
Centrale SDEGA
vista posteriore

AGSS
Station Front view
Centrale SDEGA
vista frontale

AGSS
Remote control
Controllo remoto
SDEGA

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This integrated system is used to evacuate the anaesthetic gas breathed out by the patients in the operating rooms or postoperating awakening rooms.

This system is composed of :

- electric board (local control panel);
- two vacuum pumps with lateral channels
- with double parallel flow
- one pressure transducer
- vacuum regulation valve;
- remote control system module;
- two aerating filters;

The vacuum level is adjusted by the proper regulation valve and maintained stabilized by the intermittent works of the pumps.
The transducer read and shows continuously the generated vacuum level.

Some of the advantages of this system are:

- maximum silence
- minimum overall dimensions
- safety and reliability

REMOTE CONTROL MODULE

This module, properly connected to the AGSS delta P system and installed in proximity of the operating rooms, allows to start/stop the AGSS system and control that everything is working correctly by the specific status and alarm leds. It is possible to connect a supplementary vacuum switch to monitor locally the system, for example after switch-off valves or vacuum regulators.

It is possible to install in parallel scheme more modules, the AGSS system will work if at least one of the switch is in ON position. The silencing pushbutton allows to cut off the ringer in case of alarm, if after 15 minutes the alarm conditions are still present the ringer restart automatically.

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Anesthetic gas scavenging disposal system pumps (agss/wags)

Technical Data

DESCRIPTION: The anesthetic gas scavenging disposal system pumps (figure 1) are used to evacuate anesthetic gases breathe out by patients in operating or recovery.

It is composed by:

- Electrical panel (local control panel);
- Two vacuum pumps with lateral channels and parallel streams;
- One transducer;
- One vacuum regulation valve;
- One remote controller;
- Two aeration filters.

GENERAL DATA**MAX VOLUMETRIC FLOW-RATE:**

- 50 Nm³/h
- 90 Nm³/h
- 180 Nm³/h

INLET CONNECTION: G 1" 1/2 M

OUTLET CONNECTION: G 1" 1/2 M

OPERATING TEMPERATURE: -5 ÷ +40 °C

DIMENSIONS (height x length x width) [cm]:
120 x 100 x 71

STANDARD: ISO 7396-2

THE ELECTRICAL PANEL IS COMPOSED OF:

- Main switch on/off;
- Activation switch (Pump1 / OFF / Pump2);
- Manual/Automatic modality switch;
- Pump status module;
- Control module.

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OPERATING TEMPERATURE: -5 ÷ +40 °C

STORAGE TEMPERATURE: -25 ÷ +55 °C

EQUIPPED WITH REMOTE CONTRL MODULE**Product Codes**

Anesthetic gas scavenging disposal system pumps (agss/wags)

8650-2040 ANESTHETIC GAS SCAVENGING SYSTEM 180 Nm³/h 3 Kw




8651-2040 ANESTHETIC GAS SCAVENGING SYSTEM 90 Nm³/h 1,5 Kw

8711-2040 ANESTHETIC GAS SCAVENGING SYSTEM 50 Nm³/h, 1 Kw

4489-2025M REMOTE CONTROL MODULE TO START/STOP THE SYSTEM WITH ALARM

Anesthetic gas scavenging disposal system pumps (agss/wags)

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