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### Key features

- 1 High efficiency and reliability rotary screw air compressors
- 2 Air refrigeration dryers to ensure air pressure dew point complying to **ISO 8573-1** :2010 class 4 humidity
- 3 Wide range of Air storage tank equipped with internal VITROFLEX cladding (for use with medicinal air)
- 4 Adsorption dryers for medicinal use. They allow to dehydrate compressed air, reaching a deeper dew-point than by using a refrigerated air dryer (from -20° C down to -60/-70° C). They are also equipped with filters for oil and dust removal
- 5 The **AIRTOUCH\_S** touch screen (3,5" or 5,7") control panel with PLC is designed to manage the whole air station in configurations: **2 COMPRESSORS + CYLINDERS, 3 COMPRESSORS**

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### General Description

The medical / surgical air station **AIRTOUCH\_S** is intended to supply breathable air for medical use in accordance to requirements of European Pharmacopoeia Monograph

The delta P air station has safety characteristics in order to ensure the continuity of supply also in the event of critical conditions, and anyway in normal and single fault conditions in compliance to **UNI EN ISO 7396-1** European regulatory standard

#### The main configurations managed from the control panel are:

- n°2 sources of supply (Primary and Secondary Source) composed by compressors + n.1 reserve source composed by cylinders
  - n°3 sources of supply (Primary, Secondary and Reserve Source)
- The air station is able to supply the whole design flow rate even with 2 of the 3 sources of supply out of service.

All the components of the Management Control System is integrated in a metallic box except for the electro-pneumatic block, that will be installed downstream the air vessels.

The Management and Control System monitors and displays also alarms coming from cylinders source of supply.

#### The **AIRTOUCH\_S** control panel is composed by:

Programmable Logic Controller (PLC) with colour Touch screen control panel 3.5" / 5.7"  
Modbus RTU Master communication  
Power supply 110 - 230 Vac

## MEDICAL / SURGICAL AIR STATIONS **ISO 7396-1**

### **AIRTOUCH<sub>S</sub> CONTROL PANEL MAIN FEATURES MANAGEMENT AND CONTROL SYSTEM OPERATING LOGIC**

The start/stop of compressors are controlled by the pressure transducer.  
The operating parameters of the system are programmable and the access is protected by password.  
The following operating modes are provided:

#### **AUTOMATIC MODE**

In this operating mode the PLC manages the changeover between the compressors, when one or more alarm conditions occur on the running line. Two kinds of changeover are possible: ALTERNANCE CHANGEOVER or TIMER CHANGEOVER.

#### **MAINTENANCE MODE**

In this operating mode the system by-pass all the alarm systems related to maintenance operations.

#### **MANUAL MODE**

In this operating mode the operator select the working compressor.

#### **EMERGENCY MODE**

In case of anomaly conditions (e.g. if pressure falls down under an alarm setpoint value), the system give the command of START to both the lines contemporary. Or if there is a generic anomaly on the running line, the other line is started from the system.

The control panel allows to manage the following **inputs /outputs**:

| ANALOG INPUT                                             | DIGITAL OUTPUT                                                                                         | OUTPUT                                         |
|----------------------------------------------------------|--------------------------------------------------------------------------------------------------------|------------------------------------------------|
| Pressure transducer<br>0-16 bar                          | Emergency pressure switch<br>(PE)                                                                      | Cumulative alarms<br>(operative and emergency) |
| Pressure transducer<br>0-200 bar<br>(cylinders pressure) | Compressors anomaly<br>Running of compressors<br>Compressors temperature                               | Compressors enabled                            |
| Detector CO 0-50ppm<br>[optional]                        | Filtering unit alarms<br>(Dew point, Power On,<br>DPS filter clogged)                                  |                                                |
|                                                          | Line max pressure<br>Line min pressure<br>Reserve running<br>Low pressure cylinders<br>pressure switch |                                                |

All the alarms are normally closed with alarm when contact is open.

Electro-pneumatic block

The electro-pneumatic block is composed by all the devices necessary to perform and control the process of production of medical air, both in normal and emergency conditions-

**ALARM AND MONITORING SYSTEMS**

The following alarm/warning indicators are provided:

**OPERATIVE ALARMS:**

Operative alarms are intended to indicate to technical staff that a prompt, but not immediate intervention is needed.

The following operative alarms are present:

- compressor alarm
- compressor temperature alarm
- hygrometer alarm
- differential pressure gauge alarm
- CO alarm (optional)

These indicators are provided both on the touch screen control panel and on light indicator and audible signal.

**EMERGENCY OPERATING ALARMS**

The operating emergency alarms are intended to signal all the conditions and anomalies which require an immediate response from the technical staff, such for example:

- Low or high pressure in the pipeline system
- Intervention of the emergency pressure switch (PE)
- CO alarm (optional)

**WARNING SIGNALS**

On the touch screen control panel, the following warning signals are displayed:

- Maintenance deadline warning for COMPRESSORS and FILTERING UNITS
- Warning of main power failure, at the restart of system
- Indication of switch from AUTOMATIC to MANUAL MODE

**INFORMATION SIGNALS**

- Compressor running
- Each compressor is provided with 2 hour meters. The first one is resettable from customer (by means of password), the second one is resettable only from Delta p
- Filtering unit running
- Pressure request from primary line
- Secondary line running
- Reserve source empty
- System status:
  - **No alarm:** STEADY WHITE LIGHT
  - **Operative alarm:** YELLOW FLASHING LIGHT
  - **Emergency operating alarm:** RED FLASHING LIGHT

**AIR STATION - MAIN TYPICAL LAYOUTS**

Depending on the air station configuration, the **AIRTOUCH\_S** control panel shows the following characteristics.

**2 COMPRESSORS + CYLINDERS 3rd SOURCE CONFIGURATION**

In this configuration, the **AIRTOUCH\_S** control panel is equipped to manage a system composed by:

- n.2 lines of production of medical air with compressors
  - n.1 air supply unit with cylinders
- and it works according to the above mentioned logic

**AIR STATION WITH 3 COMPRESSORS**

In this configuration the system is equipped with:

- n.2 lines of production of medical air with compressors
- n.1 additional compressor (reserve source) in parallel to the 2 primary and secondary lines with compressors

In order to be used with this configuration, in addition to the **AIRTOUCH\_S** control panel a further control panel for the reserve source (3rd compressor) must be used. When the line pressure falls down under the preset value of the pressure switch, the 3rd source control panel activates the 3rd compressor ensuring the continuity of supply to patient. When the 3rd source control panel is running, a signal is present on the **AIRTOUCH\_S** control panel. So Primary and Secondary sources are controlled by the **AIRTOUCH\_S** control panel and the 3rd Source is activated (automatically) in case of malfunctioning of Primary and Secondary Sources.

**Technical Data**

**AIR STATION GENERAL DATA**

**FLOW RATE:** From 15 to 400 Nm<sup>3</sup>/h of produced air  
**Further flow rates available on request**

**AIR OUTLET NOMINAL PRESSURE** 4 bar for medical air  
8 bar for surgical air

**COMPRESSORS WORKING PRESSURES** Systems with 10 bar and 13 bar rotary screw compressors are available as standard supply

**Further working pressures are available on request**

**AIR TANKS** 11,5 or 16 bar (depending on the compressors) marked in accordance to PED 97/23/CE standard and with an internal treatment in VITROFLEX for use with medical air.

Wide range of volumes, depending on the size of the supply system, are available

**AIR REFRIGERATION DRYERS** They ensure an air pressure dew point complying to ISO 8573-1 :2010 class 4 humidity

**FILTERING UNITS** Heatless adsorption are used, in order to reach the correct dew-point to ensure a good dehydration, from -20° C down to -60/-70° C.

**QUALITY OF PRODUCED AIR** Medical Air in accordance with requirements of European Pharmacopoeia Monograph **(SEE TABLE 1)**

**AIRTOUCH\_S CONTROL PANEL TECHNICAL DATA**

**Case dimensions (H x W x D):** 600x400x265 mm

**Case material:** Varnished metal sheet

**Working temperature:** +10°C ÷ +40 °C

**Feeding tension:** 110 - 230 Vac / 1 ph

**Frequency:** 50/60 Hz

**Analog inputs** n.3 analog inputs 4..20 mA


**Digital inputs** n.9 relay contact NC

**Outputs** n.2 free volt contact NC

**Table 1**

**MEDICAL AIR CONTAMINANTS TABLE - IN ACCORDANCE TO EUROPEAN PHARMACOPOEIA MONOGRAPH**



| CONTAMINANT                        | MAX VALUE ADMITTED              |
|------------------------------------|---------------------------------|
| OXYGEN CONCENTRATION               | ≥ 20,4% (V/V) and ≤ 21,4% (V/V) |
| TOTAL OIL CONCENTRATION            | ≤ 0,1 mg/m <sup>3</sup>         |
| CO CONCENTRATION                   | ≤ 5 ppm                         |
| CO <sub>2</sub> CONCENTRATION      | ≤ 500 ppm                       |
| WATER VAPOUR CONTENT               | ≤ 67 ppm                        |
| SO <sub>2</sub> CONCENTRATION      | ≤ 1 ppm                         |
| NO + NO <sub>2</sub> CONCENTRATION | ≤ 2 ppm                         |

MEDICAL / SURGICAL AIR STATIONS **ISO 7396-1**



Page 4 of 6

**2 COMPRESSORS + CYLINDERS 3rd SOURCE - TYPICAL LAYOUT 3D**


| LAYOUT LEGEND |                      |
|---------------|----------------------|
| Comp.         | COMPRESSOR           |
| F.U.          | FILTER UNIT          |
| DRY           | DRYER                |
| TP            | TRASDUCER (0-16 BAR) |
| PS            | PRESSURE SWITCH      |
| CYL           | 3rd SOURCE CYLINDERS |

**3 COMPRESSORS - TYPICAL LAYOUT 3D**


| LAYOUT LEGEND |                      |
|---------------|----------------------|
| Comp.         | COMPRESSOR           |
| F.U.          | FILTER UNIT          |
| DRY           | DRYER                |
| TP            | TRASDUCER (0-16 BAR) |
| PS            | PRESSURE SWITCH      |
| CYL           | 3rd SOURCE CYLINDERS |

EXAMPLE OF TYPICAL LAYOUT OF A 2 COMPRESSORS COMBINED SYSTEM MEDICAL AIR / SURGICAL AIR + 3rd SOURCE WITH CYLINDERS

| N | DESCRIPTION                              | N  | DESCRIPTION                   | N  | DESCRIPTION                                        | N  | DESCRIPTION                         |
|---|------------------------------------------|----|-------------------------------|----|----------------------------------------------------|----|-------------------------------------|
| 1 | Air Compressor                           | 6  | Pressure Gauge                | 11 | Second stage panel with by pass pressure regulator | 15 | NIST AIR Probe                      |
| 2 | Ball valve                               | 7  | Air Tank                      | 12 | Control Unit preset 3,6 - 5,7 Bar                  | 16 | Control Unit Preset 5,6 - 9,6 Bar   |
| 3 | Refrigerator Air Dryer                   | 8  | Outlet AIR 8 bar (test point) | 13 | Condensate Separator Water-Oil                     | 17 | AIRTOUCH <sub>S</sub> Control Panel |
| 4 | Filtering Unit (Special for Medical Use) | 9  | Pressure transduce (4..20 mA) | 14 | NIST Unit for AIR                                  | 18 | Cylinders 3rd source of supply      |
| 5 | Pressure relief valve                    | 10 | Emergency pressure switch     |    |                                                    |    |                                     |

EXAMPLE OF TYPICAL LAYOUT OF A 3 COMPRESSORS COMBINED SYSTEM MEDICAL AIR / SURGICAL AIR

| N | DESCRIPTION                              | N  | DESCRIPTION                   | N  | DESCRIPTION                                        | N  | DESCRIPTION                         |
|---|------------------------------------------|----|-------------------------------|----|----------------------------------------------------|----|-------------------------------------|
| 1 | Air Compressor                           | 6  | Pressure Gauge                | 11 | Second stage panel with by pass pressure regulator | 15 | NIST AIR Probe                      |
| 2 | Ball valve                               | 7  | Air Tank                      | 12 | Control Unit preset 3,6 - 5,7 Bar                  | 16 | Control Unit Preset 5,6 - 9,6 Bar   |
| 3 | Refrigerator Air Dryer                   | 8  | Outlet AIR 8 bar (test point) | 13 | Condensate Separator Water-Oil                     | 17 | AIRTOUCH <sub>S</sub> Control Panel |
| 4 | Filtering Unit (Special for Medical Use) | 9  | Pressure transduce (4..20 mA) | 14 | NIST Unit for AIR                                  | 18 | 3rd Compressor Control Panel        |
| 5 | Pressure relief valve                    | 10 | Emergency pressure switch     |    |                                                    |    |                                     |

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